

gloveon

Paloma



Nitrile Exam Gloves Powder Free, Standard Cuff

Like the purity of a white Dove, **Paloma**, our white nitrile gloves offers the elegance of touch and feel with strength and dexterity needed for the delicate procedure and jobs.



| GloveOn Paloma | | |
|---|-------------------------|-------------------------|
| Length (mm) | ≥ 230 | |
| Thickness Measurements (mm) | | |
| Palm (centre of Palm) | 0.07 ± 0.02 | |
| Finger ($13\text{mm} \pm 3\text{mm}$ from tip) | 0.09 ± 0.02 | |
| Physical Properties | Before Ageing | After Ageing |
| Tensile Strength (MPa) | ≥ 18 | ≥ 16 |
| Elongation (%) | ≥ 500 | ≥ 400 |
| Inspection Levels & AQL | Inspection Level | AQL |
| Watertightness | G1 | 1.50 |
| Physical Dimensions | S2 | 4.00 |
| Tensile Strength | S2 | 4.00 |
| Visual Inspection (Major) | S4 | 2.50 |
| Visual Inspection (Minor) | S4 | 4.00 |
| Particulate Residue | N = 5 | $\leq 2\text{mg/glove}$ |

REORDER CODE

| | |
|---------|---------|
| NTR51XS | X-SMALL |
| NTR51SS | SMALL |
| NTR51MM | MEDIUM |
| NTR51LL | LARGE |
| NTR51XL | X-LARGE |

FEATURES

- Fingertip textured
- Powder free
- Not made with natural rubber latex
- Chemo drugs tested
- Lab chemical tested
- Ambidextrous
- Standard cuff
- White colour

PACKAGING

200 gloves per box for XS to L
180 gloves per box for XL
10 boxes per carton

REGULATORY COMPLIANCE

TGA - ARTG 164563, FDA 510(k), MDD 93/42/EEC, REACH, ROHS Directive 2002/95/EC, EC 10/2011, REGULATION (EC) No 1935/2004, PPE 89/686/EEC

STANDARDS

ASTM D6319, EN 455 part 1, 2, 3 & 4, EN 1186, EN 13130, CEN/TS 14234, EN 420, EN 374 part 1, 2 & 3

MANUFACTURING ACCREDITATIONS

ISO 9001:2008
ISO 13485:2003
EN ISO 13485:2003

Chemotherapy Drugs and Concentration
(Tested for Resistance to Permeation by Chemotherapy
Drugs as per ASTM D6978-05 Test Report PN 116668)

Minimum Breakthrough
Detection Time (minutes)

| | |
|--|--------------|
| Carmustine (BCNU), 3.3mg/ml (3,300 ppm) | 15.1 minutes |
| Cisplatin, 1.0mg/ml (1,000 ppm) | >240 minutes |
| Cyclophosphamide (Cytoxan), 20.0mg/ml (20,000 ppm) | >240 minutes |
| Dacarbazine (DTIC), 10.0mg/ml (10,000 ppm) | >240 minutes |
| Doxorubicin Hydrochloride, 2.0mg/ml (2,000 ppm) | >240 minutes |
| Etoposide (Toposar), 20.0mg/ml (20,000 ppm) | >240 minutes |
| Fluorouracil, 50.0mg/ml (50,000 ppm) | >240 minutes |
| Methotrexate, 25.0mg/ml (25,000 ppm) | >240 minutes |
| Mitomycin C, 0.5mg/ml (500 ppm) | >240 minutes |
| Paclitaxel (Taxol), 6.0mg/ml (6,000 ppm) | >240 minutes |
| Thiotepa, 10.0mg/ml (10,000 ppm) | 15.4 minutes |
| Vincristine Sulfate, 1.0mg/ml (1,000 ppm) | >240 minutes |

WARNING: Carmustine and Thiotepa, at the tested concentration, degraded Paloma nitrile glove at 15.1 minutes and 15.4 minutes, respectively. The safe use of gloves in chemotherapy treatment is solely the decision of clinicians authorised to make such decision.

| Chemotherapy Drugs and Concentration (Tested for Resistance to Permeation by Chemotherapy Drugs as per ASTM D6978-05 Test Report PN 116668) | Minimum Breakthrough Detection Time (minutes) |
|---|--|
| Carmustine (BCNU), 3.3mg/ml (3,300 ppm) | Not Recommended |
| Cisplatin, 1.0mg/ml (1,000 ppm) | >240 minutes |
| Cyclophosphamide (Cytoxan), 20.0mg/ml (20,000 ppm) | >240 minutes |
| Dacarbazine (DTIC), 10.0mg/ml (10,000 ppm) | >240 minutes |
| Doxorubicin Hydrochloride, 2.0mg/ml (2,000 ppm) | >240 minutes |
| Etoposide (Toposar), 20.00mg/ml (20,000 ppm) | >240 minutes |
| Fluorouracil, 50.0mg/ml (50,000 ppm) | >240 minutes |
| Methotrexate, 25.0mg/ml (25,000 ppm) | >240 minutes |
| Mitomycin C, 0.5mg/ml (500 ppm) | >240 minutes |
| Paclitaxel (Taxol), 6.0mg/ml (6,000 ppm) | >240 minutes |
| Thiotepa, 10.0mg/ml (10,000 ppm) | Not Recommended |
| Vincristine Sulfate, 1.0mg/ml (1,000 ppm) | >240 minutes |

Warning: Don't use with Carmustine and Thiotepa

Features

- Fingertip textured
- Powder free
- Not made with natural rubber latex
- Chemo drugs tested
- Lab chemical tested
- Ambidextrous
- Standard cuff
- White colour

Sizes

- X-Small
- Small
- Medium
- Large
- X-Large

Manufacturing Accreditations

- ISO 9001:2008
- ISO 13485:2003
- EN ISO 13485:2003



MDSS GmbH
Schiffgraben 41,30175 Hannover, Germany

Certificates FDA



April 25, 2020

Hartalega NGC SDN, BHD.
Nurul Kong
Senior Manager- Quality Assurance
Kawasan Perindustrian Tanjung
Sepang, Selangor 43900
Malaysia

Re: K200581

Trade/Device Name: Paloma Nitrile Powder Free Examination Gloves (Blue)
Regulation Number: 21 CFR 880.6250
Regulation Name: Paloma Powder Free Examination Glove
Regulatory Class: Class I, reserved
Product Code: LZA, LZC, QDO
Dated: February 27, 2020
Received: March 5, 2020

Dear Nurul Kong:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the [Federal Register](#).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).


Sincerely,
Elizabeth F.
Claverie -S



CAPT Elizabeth Claverie, M.S.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Certificates CE



Certificate of CE-Registration

This is to certify that, in accordance with the Medical Device Directive 93/42/EEC, Medical Device Safety Service GmbH (MDSS) agrees to perform all duties and responsibilities as the Authorized Representative for:

**Hartalega NGC Sdn. Bhd.
No. 1, Persiaran Tanjung
Kawasan Perindustrian Tanjung
43900 Sepang, Selangor
MALAYSIA**

as stipulated and demanded by the aforementioned Directive. The German Competent Authority has allocated the medical devices of the Manufacturer registration numbers as foreseen in:

Annex A dated January 18, 2019

The Manufacturer has provided MDSS with the appropriate Declaration(s) of Conformity confirming that the medical devices fulfill the applicable requirements of Directive 93/42/EEC. In compliance with German law, a safety officer has been appointed for Germany.

2019-01-18


Lüdger Möller
President
MDSS GmbH



★ MDSS ★

Medical Device Safety Service ★



MDSS · Schiffgroben 41 · 30175 Hannover, Germany

Hartalega NGC Sdn. Bhd,
Khairunnisa Warsito
No. 1, Persiaran Tanjung
Kawasan Perindustrian Tanjung
43900 Sepang, Selangor
MALAYSIA

Schiffgroben 41
30175 Hannover, Germany

Tel: + 49 - 511 - 62 62 86 30
Fax: + 49 - 511 - 62 62 86 33

eMail: info@mdss.com
Internet: www.mdss.com

2019.01.18

Confirmation of CE Registration

Dear Khairunnisa,

It is our pleasure to enclose the new Certificate of CE-Registration for your product.

Please note that registration was performed under § 25 MPG (Medizinproduktegesetz). This is the *Federal Republic of Germany's national interpretation of Medical Device Directive 93/42/EEC*. Registration is therefore in accordance with EU legislation. We remind you that all products must meet the applicable provision of the European and national regulation before they may be placed on the market.

We are looking forward to continuing our good business relationship and wish you a successful product launch in Europe.

Best regards,

Juan Monferrer Tena
Administrative Assistant
Medical Device Safety Service GmbH

Encl.

1 Certificate of CE-Registration
1 Annex A

MDSS · Medical Device Safety Service GmbH
Handelsregister Hannover HRB 57318 · USt-IdNr. DE 177346163 · Geschäftsführen: Ludger Möller
Bankverbindungen
Sparkasse Hannover
S.W.I.F.T.: SPHDE2H
IBAN: DE24 2505 0180 0910 0792 77

Commerzbank AG, Hannover
S.W.I.F.T.: COBDEFF 350
IBAN: DE67 2504 0066 0338 8816 00



**Annex A dated January 18, 2019
 Manufacturer: Hartalega NGC Sdn. Bhd.**

| UMDNS Code Description Notified Medical Device Product Name & Catalogue Number | UMDNS Code | Risk Class | Cat. Code | Registration Number | NB No. / NB Certificate No. | NB Cert. Valid Until YYYY-MM-DD |
|--|---------------|---------------|--------------|-------------------------|--------------------------------|---------------------------------------|
| Gloves, Examination/Treatment | 11-862 | I | 10 | DE/CA09/0170/H13/001-01 | N.A. | N.A. |
| Latex Examination Gloves; Nitrile Examination Gloves; Antimicrobial Nitrile Powder Free Examination Gloves Latex Powder Free Examination Gloves Latex Powdered Examination Gloves Nitrile Powder Free Examination Gloves Antimicrobial Nitrile Powder Free Examination Gloves | | | | | | |



Certificates Testing

Testing. Development. Problem Solving.



January 5, 2019

• TEST REPORT •

PN 83672A - Amended

CHEMICAL ANALYTICAL SERVICES

Prepared For:
Hartalega SDN. BDH
Ms. Nurul Aisyah Kong
No. 7 Kawasan Perusahaan Suria
Bestari Jaya
Selangor, 45600
Malaysia

Prepared By:


Tiffany L Heller
Chemical Technician

Approved By:


Ana C. Barbur, M.S.
Manager, Chemical & Pharmaceutical Services

An A2LA Accredited Testing Laboratory — Certificate Numbers 255.01 & 255.02
ISO 9001:2000 Registered
Member of ACIL: The American Council of Independent Laboratories



ISO 9001:2000
Registered



Letters and reports are for the exclusive use of the clients to whom they are addressed and shall not be reproduced, except in full, without the written permission of Akron Rubber Development Laboratory, Inc. (ARDL). The information contained herein applies to the specific material, products or processes tested or evaluated. No warranty of any kind is herein construed or implied. The liability of ARDL, Inc. shall be limited to the amount of consideration paid for services. ARDL, Inc. is accredited by A2LA for the test methods listed on the attached scope.

www.ardl.com

2887 Gilchrist Rd. | Akron, Ohio 44305 | answers@ardl.com
Toll Free (800) 830-ARDL | Worldwide (330) 794-6600 | Fax (330) 794-6610



Testing. Development. Problem Solving.

January 5, 2019

Ms. Nurul Aisyah Kong
Hartalega SDN. BHD

Page 1 of 3 – PN 83672A - Amended

SUBJECT: Permeation testing per ASTM D 6978-05 on sample submitted by the above company. Wire Transfer.

RECEIVED: Glove sample identified as Nitrile Powder Free Examination Gloves (Blue) Code: ABLU.

TESTING CHEMOTHERAPY DRUGS:

Table 1. List of the Testing Chemotherapy Drugs, Sources, and Expiration Dates

| TESTING CHEMOTHERAPY DRUGS | DRUG SOURCE |
|----------------------------|---|
| Carmustine (BCNU) | Sigma; Lot# 038K4008; Expiration 12/2009 |
| Cisplatin | Sigma; Lot# 59H3657; Expiration 09/2009 |
| Cyclophosphamide (Cytoxan) | Sigma; Lot# 068K1131; Expiration 1/2010 |
| Dacarbazine (DTIC) | Hospira; Lot# U022223AA; Expiration 06/2010 |
| Doxorubicin Hydrochloride | Teva; Lot#07N625; Expiration 10/2009 |
| Etoposide (Toposar) | Teva; Lot# 31303976B; Expiration 9/2011 |
| Fluorouracil | APP; Lot# 203867; Expiration 03/2010 |
| Mitomycin C | Sigma; Lot# 048K1086; Expiration 01/2010 |
| Methotrexate | Hospira; Lot# U024457AA; Expiration 05/2010 |
| Paclitaxel (Taxol) | Dabur Oncology; Lot# PA08H00701; Exp. 05/2010 |
| Thiotepa | Sigma; Lot#078K1526; Expiration 12/2009 |
| Vincristine Sulfate | Hospira; Lot# U037139AA; Expiration 12/2009 |

COLLECTION MEDIA:

The collection media, which were selected, are listed in Table 2.

Table 2. Collection Media for Testing Chemotherapy Drugs

| TEST DRUG AND CONCENTRATION | COLLECTION MEDIUM |
|---|-----------------------------------|
| Carmustine (BCNU), 3.3 mg/ml (3,300 ppm) | 10% Ethanol Aqueous Solution |
| Cisplatin, 1.0 mg/ml (1,000 ppm) | Distilled Water |
| Cyclophosphamide (Cytoxan), 20 mg/ml (20,000 ppm) | Distilled Water |
| Dacarbazine (DTIC), 10.0 mg/ml (10,000 ppm) | Distilled Water |
| Doxorubicin Hydrochloride, 2.0 mg/ml (2,000 ppm) | Distilled Water |
| Etoposide (Toposar), 20.0 mg/ml (20,000 ppm) | Distilled Water |
| Fluorouracil, 50.0 mg/ml (50,000 ppm) | 9.20 pH Sodium Hydroxide Solution |
| Methotrexate, 25 mg/ml (25,000 ppm) | Distilled Water |
| Mitomycin C, 0.5 mg/ml (500 ppm) | Distilled Water |
| Paclitaxel (Taxol), 6.0 mg/ml (6,000 ppm) | 30% Methanol Aqueous Solution |
| Thiotepa, 10.0 mg/ml (10,000 ppm) | Distilled Water |
| Vincristine Sulfate, 1.0 mg/ml (1,000 ppm) | Distilled Water |

www.ardl.com

2887 Gilchrist Rd. | Akron, Ohio 44305 | answers@ardl.com
Toll Free (800) 830-ARDL | Worldwide (330) 794-6600 | Fax (330) 794-6610

Ms. Nurul Aisyah Kong
Hartalega SDN. BHD

Page 2 of 3 – PN 83672A - Amended

TESTING CONDITIONS:

| | |
|--------------------------------------|--|
| Standard Test Method Used: | ASTM D 6978-05 |
| Deviation From Standard Test Method: | Used 1" Permeation Cell |
| Analytical Method: | UV/VIS Spectrometry |
| Testing Temperature: | 35.0°C ± 2.0 |
| Collection System: | Closed Loop |
| Specimen Area Exposed: | 5.067 cm ² |
| Selected Data Points: | 25/test |
| Number of Specimens Tested: | 3/test |
| Location Sampled From: | Cuff area |
| Comments/Other Conditions: | Magnetic stir bar was used in the sampling chamber |

DETECTION METHOD OF CHEMICAL PERMEATION; UV/VIS ABSORPTION SPECTROMETRY:

Instrument: Perkin Elmer UV/VIS Spectrometer Lambda 25

UV/VIS Absorption Spectrometry was used to measure the absorbance of test chemicals, which permeated through the specimens into the collection medium. The collection medium was circulated in a closed loop at 11 ml/minute of flow rate through the testing period. Data collection was performed according to the programmed schedule by means of UV Winlab software from the Perkin Elmer Corporation. The list of the characteristic wavelengths is shown below.

Table 3. Characteristic Wavelengths used in UV/VIS Absorption Spectrometry

| TESTING CHEMOTHERAPY DRUGS | WAVELENGTH (nm) |
|---|-----------------|
| Carmustine (BCNU), 3.3 mg/ml (3,300 ppm) | 229 |
| Cisplatin, 1.0 mg/ml (1,000 ppm) | 199 |
| Cyclophosphamide (Cytoxan), 20 mg/ml (20,000 ppm) | 200 |
| Dacarbazine (DTIC), 10.0 mg/ml (10,000 ppm) | 320 |
| Doxorubicin Hydrochloride, 2.0 mg/ml (2,000 ppm) | 232 |
| Etoposide (Toposar), 20.0 mg/ml (20,000 ppm) | 205 |
| Fluorouracil, 50.0 mg/ml (50,000 ppm) | 269 |
| Methotrexate, 25 mg/ml (25,000 ppm) | 303 |
| Mitomycin C, 0.5 mg/ml (500 ppm) | 217 |
| Paclitaxel (Taxol), 6.0 mg/ml (6,000 ppm) | 231 |
| Thiotepa, 10.0 mg/ml (10,000 ppm) | 199 |
| Vincristine Sulfate, 1.0 mg/ml (1,000 ppm) | 220 |

SAMPLE CHARACTERISTICS:

Table 4. Thickness characteristics for the tested specimens: Nitrile Powder Free Examination Gloves (Blue) Code: ABLU.

| Testing Chemotherapy Drugs | Thickness (mm) | | | Average (mm) | Weight/Unit Area (g/m ²) |
|----------------------------|----------------|----------|----------|--------------|--------------------------------------|
| | Sample 1 | Sample 2 | Sample 3 | | |
| Carmustine (BCNU) | 0.080 | 0.084 | 0.082 | 0.082 | 79.2 |
| Cisplatin | 0.086 | 0.082 | 0.088 | 0.085 | 79.2 |
| Cyclophosphamide (Cytoxan) | 0.087 | 0.083 | 0.081 | 0.084 | 79.2 |
| Dacarbazine (DTIC) | 0.086 | 0.080 | 0.084 | 0.083 | 79.2 |
| Doxorubicin Hydrochloride | 0.087 | 0.084 | 0.082 | 0.084 | 79.2 |
| Etoposide (Toposar) | 0.085 | 0.088 | 0.090 | 0.088 | 79.2 |
| Fluorouracil | 0.082 | 0.082 | 0.094 | 0.086 | 79.2 |
| Methotrexate | 0.081 | 0.084 | 0.087 | 0.084 | 79.2 |
| Mitomycin C | 0.082 | 0.082 | 0.084 | 0.083 | 79.2 |
| Paclitaxel (Taxol) | 0.086 | 0.081 | 0.082 | 0.083 | 79.2 |
| Thiotepa | 0.084 | 0.083 | 0.087 | 0.085 | 79.2 |
| Vincristine Sulfate | 0.088 | 0.084 | 0.082 | 0.085 | 79.2 |

Ms. Nurul Aisyah Kong
Hartalega SDN. BHD

Page 3 of 3 – PN 83672A - Amended


RESULTS:

Table 5. Permeation Test Results on: Nitrile Powder Free Examination Gloves (Blue) Code: ABLU.

| TEST CHEMOTHERAPY DRUG AND CONCENTRATION | MINIMUM BREAKTHROUGH DETECTION TIME (Specimen 1/2/3) (Minutes) | STEADY STATE PERM. RATE (Specimen 1/2/3) ($\mu\text{g}/\text{cm}^2/\text{minute}$) | OTHER OBSERVATIONS |
|---|--|--|--------------------------------------|
| Carmustine (BCNU), 3.3 mg/ml (3,300 ppm) | 35.80 (45.31,30.46,31.63) | 1.24 (1.36,1.16,1.20) | Moderate swelling and no degradation |
| Cisplatin, 1.0 mg/ml (1,000 ppm) | No breakthrough up to 240 min. | N/A | Slight swelling and no degradation |
| Cyclophosphamide (Cytosan), 20 mg/ml (20,000 ppm) | No breakthrough up to 240 min. | N/A | Slight swelling and no degradation |
| Dacarbazine (DTIC), 10.0 mg/ml (10,000 ppm) | No breakthrough up to 240 min. | N/A | Slight swelling and no degradation |
| Doxorubicin Hydrochloride, 2.0 mg/ml (2,000 ppm) | No breakthrough up to 240 min. | N/A | Slight swelling and no degradation |
| Etoposide (Toposar), 20.0 mg/ml (20,000 ppm) | No breakthrough up to 240 min. | N/A | Slight swelling and no degradation |
| Fluorouracil, 50.0 mg/ml (50,000 ppm) | No breakthrough up to 240 min. | N/A | Slight swelling and no degradation |
| Methotrexate, 25 mg/ml (25,000 ppm) | No breakthrough up to 240 min. | N/A | Slight swelling and no degradation |
| Mitomycin C, 0.5 mg/ml (500 ppm) | No breakthrough up to 240 min. | N/A | Moderate swelling and no degradation |
| Paclitaxel (Taxol), 6.0 mg/ml (6,000 ppm) | No breakthrough up to 240 min. | N/A | Moderate swelling and no degradation |
| Thiotepa, 10.0 mg/ml (10,000 ppm) | 85.48 (105.33,75.52,75.60) | 1.32 (1.33,1.32,1.32) | Slight swelling and no degradation |
| Vincristine Sulfate, 1.0 mg/ml (1,000 ppm) | No breakthrough up to 240 min. | N/A | Slight swelling and no degradation |




Tiffany L. Heller
Chemical Technician, Chemical Services
AKRON RUBBER DEVELOPMENT LABORATORY, INC.



Ana C. Barbur, M.S.
Chemical and Pharmaceutical Services

Certificates Testing

| MATERIAL SAFETY DATA SHEET | |  | | |
|---|---------------------------------|---|------------|------------|
| SECTION 1: PRODUCT IDENTIFICATION | | | | |
| NAME Hartalega Sdn. Bhd. | | ADDRESS C-G-9, Jalan Dataran SD1, Dataran SD, PJU 9, Bandar Sri Damansara, 52200 Kuala Lumpur | | |
| TELEPHONE NUMBER (603) 62771733 | | DATE PREPARED October 15, 2014 | | |
| COMMON NAME (USED ON LABEL) Nitrile Powder Free Examination Gloves | | CHEMICAL FAMILY Carboxylated Butadiene Acrylonitrile Polymer Latex | | |
| APPLICATION Medical and Dental | | TRADENAME & SYNONYM GLOVEON PALOMA (NTR51) NITRILE POWDER FREE EXAMINATION GLOVES (WHITE) | | |
| SECTION 2: HAZARDOUS INGREDIENTS | | | | |
| HAZARDOUS COMPONENT | CAS # | %(WT) | TLV | PEL |
| N/A | N/A | N/A | N/A | N/A |
| <small>PEL : Permissible Exposure Limit established by Occupational Safety and Health Administration (OSHA). TLV : Threshold Limit Value established by the American Conference of Governmental Industrial Hygienists, 1987-1988</small> | | | | |
| SECTION 3: COMPOSITION/ INFORMATION ON INGREDIENTS | | | | |
| CHEMICAL COMPOSITION All chemicals used are non-toxic/non-hazardous. Butadiene-Acrylonitrile Latex, Sodium Dodecylbenzenesulfonate, Sulphur, Titanium Dioxide, Zinc-n-dibutylidithiocarbamate, Zinc Oxide, Paraffin Wax Emulsion | | | | |
| SECTION 4: FIRST AID MEASURE | | | | |
| If reaction in the form of skin irritation is noticed, remove gloves immediately and wash affected part with saline water. If there is no relief, seek medical reactions. | | | | |
| SECTION 5: FIRE FIGHTING MEASURE | | | | |
| FLASHPOINT | AUTOIGNITION TEMPERATURE | FLAMMABLE LIMITS IN AIR | | |
| N/A | N/A | N/A | | |
| EXTINGUISHING MEDIA Chemical foam and dry chemical may be used | | | | |
| FIRE-FIGHTING PROCEDURES Use standard procedures for combustion material fires, including approved self-contained breathing apparatus | | | | |
| FIRE AND EXPLOSION HAZARDS No fire or explosion hazards are associated with these products. They will melt at elevated temperatures | | | | |
| SECTION 6: ACCIDENTAL RELEASE MEASURES | | | | |
| BIOCOMPATABILITY The chemical formulation of the gloves and surface lubrication materials does not contain any substances normally known to be harmful to the user or to any person with whom the gloves come into contact. | | | | |
| MEDICAL CONDITIONS GENERALLY AGGRAVATED BY EXPOSURE: Nitrile Powder Free Gloves are not expected to cause any adverse health effects. | | | | |
| SECTION 7: HANDLING AND STORAGE | | | | |
| PRECAUTIONS TO BE TAKEN IN HANDLING AND STORAGE Store in a dry, cool and ventilated area. Do not store above 104°F (40°C). Shield open box from direct sunlight, fluorescent lighting and x-rays. Improper storage will decrease usable life. | | | | |

| SECTION 8: EXPOSURE CONTROLS/ PERSONAL PROTECTION | | | | | |
|--|----------------------------|---|---|-------------------------------------|--------------------|
| EYE PROTECTION Not necessary under conditions of intended use. | | | SKIN PROTECTION Not necessary under conditions of intended use. | | |
| RESPIRATORY PROTECTION Not necessary under conditions of intended use. | | | VENTILATION Not necessary under conditions of intended use. | | |
| STEPS TO BE TAKEN IN CASE MATERIAL IS LEAKED OR SPILLED These products are solid articles and are not subject to leaks or spills. | | | | | |
| SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES | | | | | |
| APPEARANCE/ODOR Ambidextrous, Beaded Cuff, Micro-textured, Chlorinated, Powder Free, White. | | | | | |
| DIMENSION | EXTRA-SMALL | SMALL | MEDIUM | LARGE | EXTRA-LARGE |
| Length (mm) | Minimum 230 (same for all) | | | | |
| Width (mm) | 76±4 | 86±4 | 98±4 | 107±4 | 115±4 |
| THICKNESS (mm) – SINGLE WALL MEASUREMENT (same for all) | | | | | |
| Finger (mm) | 0.09 ± 0.02 | | | | |
| Palm (mm) | 0.07 ± 0.02 | | | | |
| TENSILE PROPERTIES | | UNAGED | | AGED | |
| Tensile Strength (Mpa) | | Min. 14.0 MPa | | Min. 14.0 MPa | |
| Ultimate Elongation (%) | | Min. 500% | | Min. 400% | |
| SECTION 10: STABILITY AND REACTIVITY | | | | | |
| BOILING POINT N/A | | VAPOR PRESSURE (mm Hg) N/A | | VAPOR DENSITY (air=1) N/A | |
| SPECIFIC GRAVITY (water=1) N/A | | SOLUBILITY IN WATER Insoluble | | % VOLATILE BY VOLUME N/A | |
| EVAPORATION RATE N/A | | | VISCOSITY N/A | | |
| SECTION 11: TOXICOLOGICAL INFORMATION | | | | | |
| STABILITY Stable | | | CONDITIONS TO AVOID Does not apply | | |
| INCOMPATIBILITY (MATERIALS TO AVOID) High polar solvent like methyl ethyl ketone, acetone. | | | | | |
| HAZARDOUS DECOMPOSITION PRODUCTS: In a fire, these products may produce a black smoke. Carbon Dioxide, Carbon Monoxide, Oxides of Nitrogen, aromatic/aliphatic hydrocarbons. | | | | | |
| HAZARDOUS POLYMERIZATION Will not occur | | | | | |
| SECTION 12: ECOLOGICAL INFORMATION | | | | | |
| NA | | | | | |
| SECTION 13: DISPOSAL CONSIDERATION | | | | | |
| WASTE DISPOSAL METHOD Consult current local, state and federal regulations for proper disposal methods. | | | | | |
| SECTION 14: TRANSPORT INFORMATION | | | | | |
| NA | | | | | |
| SECTION 15: REGULATORY INFORMATION | | | | | |
| NA | | | | | |
| SECTION 16: OTHER INFORMATION | | | | | |
| RECOMMENDED USE AND RESTRICTION The Nitrile Powder Free Gloves is a Single Use device | | | | | |