TOP GLOVE

TOP QUALITY, TOP EFFICIENCY

THE WORLD'S LARGEST MANUFACTURER OF GLOVES

NITRILE GLOVE

POWDERED & POWDER-FREE



Functional Benefits:

- Protection from unwanted or dangerous substances
- · Beaded cuff makes donning easy and helps prevent roll back
- Superior strength with better puncture resistance
- · Full textured enhances wet and dry grip
- Thinner gauger improves tactile sensitivity
- Custom design enhances comfort and fit
- Provide an alternative solution for individuals who are allergic to natural rubber latex



NITRILE GLOVE

NITRILE GLOVE

POWDERED & POWDER-FREE

Quality Standards

- Conforms to ASTM D6319 and EN455 Standards
- Manufactured under QSR (GMP), ISO 9001:2015 and ISO 13485:2016 Quality Management System

Glove Sizes

- Extra-Small, Small, Medium, Large, Extra-Large
- · Size of gloves shall be marked in the check box on the shipping carton with black ink

Product Specifications

Туре	:	Powdered & Powder-Free, Non-sterile
Material	:	100% Synthetic Nitrile Latex
Colour	:	Blue, White, Green, Pink, Light Purple, Black, Red
Design & Features	:	Powdered:
		Ambidextrous, finger textured or palm textured surface,
		beaded cuff, USP grade absorbable cornstarch
		Powder-Free:
		Polymer coated or online single chlorinated, offline double chlorinated, ambidextrous, finger textured or palm textured surface, beaded cuff
Storage	:	The gloves shall maintain their properties when stored in a dry condition. Avoid direct sunlight

Shelf-life : 5 years from the date of manufacturing

Physical Dimensions

Dimensions	Standards				
	Top Glove	ASTM D3578	EN 455		
Length (mm)	Min 230, Min 240, 300 ± 10	Min 220 (XS, S) Min 230 (M, L, XL)	Min 240		
Palm Width (mm)					
• XS	76 ± 3	70 ± 10	≤ 80		
• S	84 ± 3	80 ± 10	80 ± 10		
• M	94 ± 3	95 ± 10	95 ± 10		
۰L	105 ± 3	110 ± 10	110 ± 10		
• XL	113 ± 3	120 ± 10	≥ 110		
Thickness : Single Wall (mm)					
Finger	Min 0.05	Min 0.05	N/A		
• Palm	Min 0.05	Min 0.05	N/A		

If you have any enquiries on our products, please contact :

TOP GLOVE SALES & CORPORATE OFFICE

Level 21, Top Glove Tower, 16, Persiaran Setia Dagang, Setia Alam, Seksyen U13, 40170 Shah Alam, Selangor D.E., Malaysia. Tel : +603-3362 3098 : +603-3362 5096 (Sales) Fax : sales@topglove.com.my Email

Physical Properties

Property	ASTM D6319	EN 455
Tensile Strength (MPa)		
Before Aging	Min 14	N/A
After Aging	Min 14	N/A
Elongation at Break (%)		
Before Aging	Min 500	N/A
After Aging	Min 400	N/A
Median Force at Break (N)		
Before Aging	N/A	Min 6
After Aging	N/A	Min 6

International Quality Certificate Awarded:





https://www.facebook.com/mytopglove

Top Glove Website





f

