





Test Report No.: QDHL2004003051MD-01 Date: SEP.22,2020 Page: 1 of 5

SHANDONG INTCO MEDICAL PRODUCTS CO., LTD NO.9888, QIWANG ROAD, NAOSHAN INDUSTRY PARK, QINGZHOU, SHANDONG, CHINA

THIS REPORT IS TO SUPERSEDE TEST REPORT NO.: QDHL2004003051MD, DATE: APR.30,2020. THE ORIGINAL REPORT SHALL BECOME INVALID AS OF THE DATE OF ISSUANCE OF THIS REPORT.

The following sample(s) was/were submitted and identified by the client as:

Sample Description : SYNMAX VINYL EXAM GLOVESA

Lot No. : NOT PROVIDED

Lot Size : NOT PROVIDED

Sample Quantity : 80PCS

Sample Receiving Date : APR.17,2020

Testing Period : APR.17,2020 TO APR.30,2020

Test Performed : SELECTED TEST(S) AS REQUESTED BY APPLICANT

Test Requested : BS EN 455-2:2015 MEDICAL GLOVES FOR SINGLE USE – PART 2:

REQUIREMENTS AND TESTING FOR PHYSICAL PROPERTIES

(CLAUSE 4.2, 4.3, 5.2, 5.3)

Test Result(s) : PLEASE REFER TO THE FOLLOWING PAGE(S)

Remark: Unless otherwise stated the results shown in this test report refer only to the sample(s) tested. This document cannot be used for publicity, without prior written approval of the SGS.

SGS-CSTC Standards Technical Services (Qingdao) Co., Ltd.

Jessica Gao

Approved Signatory



QDHL2004003051MD-0



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com/en/Terms-and-Conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction lssues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

QD

7389996

www.sgsgroup.com.cn e sgs.china@sgs.com







Test Report

No.: QDHL2004003051MD-01

Date: SEP.22,2020

Page: 2 of 5

Test Conducted:

BS EN 455-2:2015 Medical gloves for single use – part 2: Requirements and testing for physical properties

Number of test sample		26 Pieces	
Type	6 9	Examination/procedure gloves c)	
Size : Examination/procedure gloves: M		Examination/procedure gloves: M	

Clause	Test Items	Result	Note
4 Dimensions		Z (23 C3)	F
4.2	Length	Pass	# 1
4.3	Width	Pass	#1
5 Strength		234	-15"
5.2 Force at break		Pass	# 2
5.3 Force at break after challenge testing		Pass	#2

Note	:	#1	See result 1
-6		#2	See result 2



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com/en/Terms-and-Conditions/Terms

QD

7389997

www.sgsgroup.com.cn e sgs.china@sgs.com







Test Report

No.: QDHL2004003051MD-01

Date: SEP.22,2020

Page: 3 of 5

Test Result:

Result 1. Length and width

No.	Size: M		
140.	Length (mm)	Width (mm)	
1	240	97	
2	241	97	
3	235	97	
4	239	96	
5	235	96	
6	236	96	
7 9	241	97	
8	237	97	
9 241		97	
10	241	96	
11 4	240	97	
12	240	98	
13	240	96	
Standard requirement	≥240	95±10	
Median value	240	97	



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

QD

7389998

www.sgsgroup.com.cn e sgs.china@sgs.com 04







Test Report No.: QDHL2004003051MD-01 Date: SEP.22,2020

Result 2. Force at break

9 40 0	Size: M			
No.	Force at break (N)			
G 9"	Before aging	After aging		
1	3.8	3.4		
2	3.6	3.6		
3	4.3	3.6		
4	3.8	3.7		
5	3.6	3.6		
6	3.8	3.9		
7	3.7	3.7		
8	4.1	3.9		
9	4.0	3.6		
10	3.6	3.9		
11	3.6	3.6		
12	3.7	3.8		
13	3.6	3.5		
Standard ≥3.6		≥3.6		
Median value	3.7	3.6		

Remark:

(1)The declaration of conformity is only based on the actual value of laboratory activity, measurement uncertainty of the results not take into account.

(2) The marked part is the modification information.



Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com/en/Terms-and-Conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.



Page: 4 of 5








Test Report

No.: QDHL2004003051MD-01

Date: SEP.22,2020

Page: 5 of 5

Sample Photo:



SGS authenticate the photo on original report only

End of Report





Unless otherwise agreed in writing, this document is issued by the Company subject to its General Conditions of Service printed overleaf, available on request or accessible at http://www.sgs.com/en/Terms-and-Conditions.aspx and, for electronic format documents, subject to Terms and Conditions for Electronic Documents at http://www.sgs.com/en/Terms-and-Conditions/Terms-e-Document.aspx. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained hereon reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents. This document cannot be reproduced except in full, without prior written approval of the Company. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law.

QD

7390000

www.sgsgroup.com.cn e sgs.china@sgs.com